

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **Product datasheet for TA506233**

## PDGF Receptor beta (PDGFRB) Mouse Monoclonal Antibody [Clone ID: OTI1D2]

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: OTI1D2

Applications: FC, IF, WB

Recommend Dilution: WB 1:2000, IF 1:100, FLOW 1:50

Reactivity: Human
Host: Mouse
Isotype: IgG2a

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human PDGFRB(NP\_002600) produced in HEK293T

cell.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 123.8 kDa

**Gene Name:** platelet derived growth factor receptor beta

Database Link: NP 002600 Entrez Gene 5159 Human

**Background:** This gene encodes a cell surface tyrosine kinase receptor for members of the platelet-derived

growth factor family. These growth factors are mitogens for cells of mesenchymal origin. The

identity of the growth factor bound to a receptor monomer determines whether the functional receptor is a homodimer or a heterodimer, composed of both platelet-derived growth factor receptor alpha and beta polypeptides. This gene is flanked on chromosome 5 by the genes for granulocyte-macrophage colony-stimulating factor and macrophage-colony

stimulating factor receptor; all three genes may be implicated in the 5-q syndrome. A translocation between chromosomes 5 and 12, that fuses this gene to that of the translocation, ETV6, leukemia gene, results in chronic myeloproliferative disorder with

eosinophilia. [provided by RefSeq, Jul 2008]

Synonyms: CD140B; IBGC4; IMF1; JTK12; PDGFR; PDGFR-1; PDGFR1






Protein Families: Druggable Genome, ES Cell Differentiation/IPS, Protein Kinase, Transmembrane


**Protein Pathways:** Calcium signaling pathway, Colorectal cancer, Cytokine-cytokine receptor interaction, Focal adhesion, Gap junction, Glioma, MAPK signaling pathway, Melanoma, Pathways in cancer,

Prostate cancer, Regulation of actin cytoskeleton


## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY PDGFRB ([RC206377], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-PDGFRB. Positive lysates [LY400921] (100ug) and [LC400921] (20ug) can be purchased separately from OriGene.




Anti-PDGFRB mouse monoclonal antibody (TA506233) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY PDGFRB ([RC206377]).



Immunofluorescent staining of HeLa cells using anti-PDGFRB mouse monoclonal antibody (TA506233).





Flow cytometric Analysis of living SF295 cells, using anti-PDGFRB antibody (TA506233), (Red), compared to a nonspecific negative control antibody, (Blue).