

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA506343

CD64 (FCGR1A) Mouse Monoclonal Antibody [Clone ID: OTI3D3]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI3D3

Applications: FC, IF, IHC, WB

Recommend Dilution: WB 1:4000, IHC 1:150, IF 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human FCGR1A(NP_000557) produced in HEK293T

cell.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 42.5 kDa

Gene Name: Fc fragment of lgG receptor la

Database Link: NP 000557 Entrez Gene 2209 Human

Background: This gene encodes a protein that plays an important role in the immune response. This

protein is a high-affinity Fc-gamma receptor. The gene is one of three related gene family

members located on chromosome 1. [provided by RefSeq, Jul 2008]

Synonyms: CD64; CD64A; FCRI; IGFR1

Protein Families: Transmembrane


Protein Pathways: Fc gamma R-mediated phagocytosis, Hematopoietic cell lineage, Systemic lupus

erythematosus







Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY FCGR1A ([RC207487], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-FCGR1A. Positive lysates [LY400193] (100ug) and [LC400193] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (10ug) from 4 different cell lines by using anti-FCGR1A monoclonal antibody at 1:200 dilution.




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human ovary tissue using anti-FCGR1A mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, TA506343)






Immunohistochemical staining of paraffinembedded Human lymphoma tissue using anti-FCGR1A mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, TA506343)



Anti-FCGR1A mouse monoclonal antibody (TA506343) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY FCGR1A ([RC207487]).



Flow cytometric analysis of living 293T cells transfected with FCGR1A overexpression plasmid ([RC207487]), Red)/empty vector ([PS100001], Blue) using anti-FCGR1A antibody (TA506343). Cells incubated with a non-specific antibody (Green) were used as isotype control (1:100).