

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA506581

Tubulin (TUBA1B) Mouse Monoclonal Antibody [Clone ID: OTI3H2]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI3H2

Applications: WB

Recommend Dilution: WB 1:2000

Reactivity: Human Host: Mouse

Isotype: IgG2a

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human TUBA1B(NP_006073) produced in HEK293T

cell.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 50 kDa

Gene Name: tubulin alpha 1b

Database Link: NP 006073 Entrez Gene 10376 Human

Synonyms: K-ALPHA-1


Protein Families: Druggable Genome

Protein Pathways: Gap junction, Pathogenic Escherichia coli infection





Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY TUBA1B ([RC205008], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TUBA1B. Positive lysates [LY401831] (100ug) and [LC401831] (20ug) can be purchased separately from OriGene.