

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA506769

RPIP8 (RUNDC3A) Mouse Monoclonal Antibody [Clone ID: OTI3C6]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI3C6
Applications: IHC, WB

Recommend Dilution: WB 1:1000, IHC: 1:150

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human RUNDC3A(NP_006686) produced in

HEK293T cell.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 45.5 kDa


Gene Name: RUN domain containing 3A

Database Link: NP 006686 Entrez Gene 10900 Human


Synonyms: RAP2IP; RPIP-8; RPIP8



Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY RUNDC3A ([RC200211], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-RUNDC3A. Positive lysates [LY402004] (100ug) and [LC402004] (20ug) can be purchased separately from OriGene.



Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human colon tissue using anti-RUNDC3A mouse monoclonal antibody. (TA506769)





Immunohistochemical staining of paraffinembedded Human pancreas tissue within the normal limits using anti-RUNDC3A mouse monoclonal antibody. (TA506769)