

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA506840

p60 CAF1 (CHAF1B) Mouse Monoclonal Antibody [Clone ID: OTI4F7]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI4F7

Applications: IF, IHC, WB

Recommend Dilution: WB 1:200~1000, IHC 1:150, IF 1:100

Reactivity: Human, Monkey, Mouse, Rat, Dog

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human CHAF1B(NP_005432) produced in HEK293T

cell

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 61.3 kDa

Gene Name: chromatin assembly factor 1 subunit B

Database Link: NP 005432 Entrez Gene 110749 MouseEntrez Gene 288242 RatEntrez Gene 487750

<u>DogEntrez Gene 695206 MonkeyEntrez Gene 8208 Human</u>

Background: Chromatin assembly factor I (CAF-I) is required for the assembly of histone octamers onto

newly-replicated DNA. CAF-I is composed of three protein subunits, p50, p60, and p150. The protein encoded by this gene corresponds to the p60 subunit and is required for chromatin assembly after replication. The encoded protein is differentially phosphorylated in a cell cycle-dependent manner. In addition, it is normally found in the nucleus except during mitosis, when it is released into the cytoplasm. This protein is a member of the WD-repeat

HIR1 family and may also be involved in DNA repair. [provided by RefSeq, Jul 2008]


Synonyms: CAF-1; CAF-IP60; CAF1; CAF1A; CAF1P60; MPHOSPH7; MPP7

Protein Families: Druggable Genome, Transcription Factors







Product images:



Anti-CHAF1B mouse monoclonal antibody (TA506840) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY CHAF1B ([RC202858]).




Immunohistochemical staining of paraffinembedded Human lymphoma tissue using anti-CHAF1B mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, TA506840)




Immunohistochemical staining of paraffinembedded Human lymph node tissue within the normal limits using anti-CHAF1B mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, TA506840)






Immunohistochemical staining of paraffinembedded Carcinoma of Human prostate tissue using anti-CHAF1B mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, TA506840)



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY CHAF1B ([RC202858], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-CHAF1B. Positive lysates [LY401669] (100ug) and [LC401669] (20ug) can be purchased separately from OriGene.



Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-CHAF1B monoclonal antibody (HepG2: human; HeLa: human; SVT2: mouse; A549: human; COS7: monkey; Jurkat: human; MDCK: canine; PC12: rat; MCF7: human).