

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA506850

IL1 alpha (IL1A) Mouse Monoclonal Antibody [Clone ID: OTI2F11]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI2F11

Applications: IF, IHC, WB

Recommend Dilution: WB 1:2000, IHC 1:150, IF 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human IL1A(NP_000566) produced in HEK293T

cell.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 30.4 kDa

Gene Name: interleukin 1 alpha

Database Link: NP 000566 Entrez Gene 3552 Human

Background: The protein encoded by this gene is a member of the interleukin 1 cytokine family. This

cytokine is a pleiotropic cytokine involved in various immune responses, inflammatory processes, and hematopoiesis. This cytokine is produced by monocytes and macrophages as a proprotein, which is proteolytically processed and released in response to cell injury, and thus induces apoptosis. This gene and eight other interleukin 1 family genes form a cytokine gene cluster on chromosome 2. It has been suggested that the polymorphism of these genes is associated with rheumatoid arthritis and Alzheimer's disease. [provided by RefSeq, Jul

2008]

Synonyms: IL-1A; IL1; IL1-ALPHA; IL1F1

Protein Families: Druggable Genome, Secreted Protein







Protein Pathways:

Apoptosis, Cytokine-cytokine receptor interaction, Graft-versus-host disease, Hematopoietic cell lineage, MAPK signaling pathway, Prion diseases, Type I diabetes mellitus


Product images:



HEK293T cells lysate (5ug, left lane) and full length human recombinant protein of human IL1A (NP_000566) produced in HEK293T cell (5ug, right lane)were separated by SDS-PAGE and immunoblotted with anti-IL1A.




Immunohistochemical staining of paraffinembedded Human pancreas tissue within the normal limits using anti-IL1A mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, TA506850)



Immunohistochemical staining of paraffinembedded Carcinoma of Human pancreas tissue using anti-IL1A mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, TA506850)





Anti-IL1A mouse monoclonal antibody (TA506850) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY IL1A ([RC202084]).