

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA506977**

## **Necdin (NDN) Mouse Monoclonal Antibody [Clone ID: OTI1B1]**

### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI1B1
Applications: WB

Reactivity: WB 1:1000
Human

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human NDN(NP\_002478) produced in HEK293T

cell

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 35.9 kDa

**Gene Name:** necdin, MAGE family member

Database Link: NP 002478 Entrez Gene 4692 Human

Background: This intronless gene is located in the Prader-Willi syndrome deletion region. It is an imprinted

gene and is expressed exclusively from the paternal allele. Studies in mouse suggest that the protein encoded by this gene may suppress growth in postmitotic neurons. [provided by

RefSeq, Jul 2008]


Synonyms: HsT16328; PWCR

**Protein Families:** Druggable Genome, Transcription Factors





# **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY NDN ([RC200448], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-NDN. Positive lysates [LY400887] (100ug) and [LC400887] (20ug) can be purchased separately from OriGene.