

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA507230

HORMAD2 Mouse Monoclonal Antibody [Clone ID: OTI2E2]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI2E2
Applications: WB

Recommend Dilution: WB 1:4000

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human HORMAD2(NP_689723) produced in

HEK293T cell.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 35.1 kDa


Gene Name: HORMA domain containing 2

Database Link: NP 689723 Entrez Gene 150280 Human

Synonyms: CT46.2



Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY HORMAD2 ([RC205453], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-HORMAD2. Positive lysates [LY407512] (100ug) and [LC407512] (20ug) can be purchased separately from OriGene.