

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA800037

USP12 Mouse Monoclonal Antibody [Clone ID: OTI2A6]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI2A6

Applications: WB

Recommend Dilution: WB 1:2000

Reactivity: Human Host: Mouse

Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 61-251 of human

USP12 (NP_872294) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 42.7 kDa

Gene Name: ubiquitin specific peptidase 12


Database Link: NP 872294 Entrez Gene 219333 Human

Synonyms: UBH1; USP12L1

Protein Families: Druggable Genome, Protease



Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY USP12 ([RC205355], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-USP12. Positive lysates [LY403634] (100ug) and [LC403634] (20ug) can be purchased separately from OriGene.