

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA800429

ACADS Mouse Monoclonal Antibody [Clone ID: OTI1D2]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI1D2
Applications: IHC, WB

Recommend Dilution: WB 1:500, IHC 1:150

Reactivity: Human, Monkey

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 34-260 of human

ACADS (NP_000008) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 41.7 kDa

Gene Name: acyl-CoA dehydrogenase, C-2 to C-3 short chain

Database Link: NP 000008 Entrez Gene 698196 MonkeyEntrez Gene 35 Human

Background: This gene encodes a a tetrameric mitochondrial flavoprotein, which is a member of the acyl-

CoA dehydrogenase family. This enzyme catalyzes the initial step of the mitochondrial fatty acid beta-oxidation pathway. Mutations in this gene have been associated with Short Chain

Acyl-CoA Dehydrogenase Deficiency. [provided by RefSeq, Jul 2008]

Synonyms: ACAD3; SCAD

Protein Families: Druggable Genome


Protein Pathways: Butanoate metabolism, Fatty acid metabolism, Metabolic pathways, Valine, leucine and

isoleucine degradation







Product images:



Immunohistochemical staining of paraffinembedded Human pancreas tissue within the normal limits using anti-ACADS mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA800429)




Immunohistochemical staining of paraffinembedded Human liver tissue within the normal limits using anti-ACADS mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA800429)



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ACADS ([RC206855], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ACADS. Positive lysates [LY400002] (100ug) and [LC400002] (20ug) can be purchased separately from OriGene.





Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-ACADS monoclonal antibody (HepG2: human; HeLa: human; SVT2: mouse; A549: human; COS7: monkey; Jurkat: human; MDCK: canine; PC12: rat; MCF7: human).