

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA800559

CCDC93 Mouse Monoclonal Antibody [Clone ID: OTI6E11]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI6E11
Applications: IHC, WB

Reactivity: WB 1:2000, IHC 1:150 **Human, Monkey, Dog**

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 293-631 of human

CCDC93 (NP_061917) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 73 kDa


Gene Name: coiled-coil domain containing 93

Database Link: NP 061917 Entrez Gene 476119 DogEntrez Gene 694513 MonkeyEntrez Gene 54520 Human


Synonyms: FLJ10996; FLJ25197; MGC13033



Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY CCDC93 ([RC206101], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-CCDC93. Positive lysates [LY412751] (100ug) and [LC412751] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-CCDC93 monoclonal antibody (HepG2: human; HeLa: human; SVT2: mouse; A549: human; COS7: monkey; Jurkat: human; MDCK: canine; PC12: rat; MCF7: human).






Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human endometrium tissue using anti-CCDC93 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA800559)



Immunohistochemical staining of paraffinembedded Human lymph node tissue within the normal limits using anti-CCDC93 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA800559)



Immunohistochemical staining of paraffinembedded Human lymphoma tissue using anti-CCDC93 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA800559)