

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA800606

ZNF307 (ZKSCAN4) Mouse Monoclonal Antibody [Clone ID: OTI6C7]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI6C7

Applications: WB

Recommend Dilution: WB 1:2000

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 1-316 of human

ZKSCAN4 (NP_061983) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 61.4 kDa

Gene Name: zinc finger with KRAB and SCAN domains 4


Database Link: NP 061983 Entrez Gene 387032 Human

Synonyms: P1P373C6; ZNF307; ZNF427; ZSCAN36

Protein Families: Transcription Factors



Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ZKSCAN4 ([RC204276], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ZKSCAN4. Positive lysates [LY402737] (100ug) and [LC402737] (20ug) can be purchased separately from OriGene.