

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA801456

Thyroglobulin (TG) Mouse Monoclonal Antibody [Clone ID: OTI3E12]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI3E12
Applications: IHC, WB

Recommend Dilution: WB 1:2000, IHC 1:150

Reactivity: Human
Host: Mouse
Isotype: IgG2b

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 2511-2768 of human

TG (NP_003226) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Gene Name: thyroglobulin

Database Link: NP 003226 Entrez Gene 7038 Human

Background: Thyroglobulin (Tg) is a glycoprotein homodimer produced predominantly by the thryroid

gland. It acts as a substrate for the synthesis of thyroxine and triiodothyronine as well as the storage of the inactive forms of thyroid hormone and iodine. Thyroglobulin is secreted from the endoplasmic reticulum to its site of iodination, and subsequent thyroxine biosynthesis, in the follicular lumen. Mutations in this gene cause thyroid dyshormonogenesis, manifested as

goiter, and are associated with moderate to severe congenital hypothyroidism.

Polymorphisms in this gene are associated with susceptibility to autoimmune thyroid diseases (AITD) such as Graves disease and Hashimoto thryoiditis. [provided by RefSeq, Nov

2009]

Synonyms: AITD3; TGN


Protein Families: Druggable Genome

Protein Pathways: Autoimmune thyroid disease






Product images:



Human recombinant protein fragment corresponding to amino acids 2511-2768 of human TG (NP_003226) produced in E.coli.



Immunohistochemical staining of paraffinembedded Human thyroid tissue within the normal limits using anti-TG mouse monoclonal antibody. (TA801456; heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris, pH8.5, 120°C for 3min)



Immunohistochemical staining of paraffinembedded Carcinoma of Human thyroid tissue using anti-TG mouse monoclonal antibody. (TA801456; heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris, pH8.5, 120°C for 3min)