

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA804501

Cardiac Troponin I (TNNC1) Mouse Monoclonal Antibody [Clone ID: OTI8F8]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI8F8
Applications: IHC, WB

Recommend Dilution: WB 1:2000, IHC 1:150

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human TNNC1 (NP_003271) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 18.2 kDa

Gene Name: troponin C1, slow skeletal and cardiac type

Database Link: NP 003271 Entrez Gene 7134 Human

Background: Troponin is a central regulatory protein of striated muscle contraction, and together with

tropomyosin, is located on the actin filament. Troponin consists of 3 subunits: Tnl, which is the inhibitor of actomyosin ATPase; TnT, which contains the binding site for tropomyosin; and TnC, the protein encoded by this gene. The binding of calcium to TnC abolishes the inhibitory action of Tnl, thus allowing the interaction of actin with myosin, the hydrolysis of

ATP, and the generation of tension. Mutations in this gene are associated with

cardiomyopathy dilated type 1Z. [provided by RefSeq, Oct 2008]

Synonyms: CMD1Z; CMH13; TN-C; TNC; TNNC


Protein Pathways: Calcium signaling pathway, Cardiac muscle contraction, Dilated cardiomyopathy,

Hypertrophic cardiomyopathy (HCM)






Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY TNNC1 ([RC207014], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TNNC1. Positive lysates [LY401131] (100ug) and [LC401131] (20ug) can be purchased separately from OriGene.



Immunohistochemical staining of paraffinembedded Human adult heart tissue using anti-TNNC1 mouse monoclonal antibody. (Heatinduced epitope retrieval by 1 mM EDTA in 10mM Tris, pH8.5, 120°C for 3min, TA804501)