

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA804928

SAV1 Mouse Monoclonal Antibody [Clone ID: OTI7C7]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI7C7

Applications: WB

Recommend Dilution: WB 1:2000

Reactivity: Human Host: Mouse

Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 167-383 of human SAV1

(NP_068590) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 44.5 kDa

Gene Name: salvador family WW domain containing protein 1

Database Link: NP 068590 Entrez Gene 60485 Human

Background: WW domain-containing proteins are found in all eukaryotes and play an important role in the

regulation of a wide variety of cellular functions such as protein degradation, transcription, and RNA splicing. This gene encodes a protein with two WW domains, a SARAH domain, and a coiled-coil region and is ubiquitously expressed in adult tissues. This protein binds to MST1 (mammalian sterile 20-like kinase 1) and promotes MST1-induced apoptosis. It has also been shown to bind to HAX1 (hematopoietic cell-specific protein 1 (HS1)-associated protein X-1) and to attenuate the anti-apoptotic effects of HAX1. Studies in human and mouse suggest


this gene acts as a tumor suppressor. [provided by RefSeq, Aug

Synonyms: SAV; WW45; WWP4





Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY SAV1 ([RC204389], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-SAV1. Positive lysates [LY402881] (100ug) and [LC402881] (20ug) can be purchased separately from OriGene.