

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA805006**

## NR0B1 Mouse Monoclonal Antibody [Clone ID: OTI5F5]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI5F5

Applications: IHC, WB

Recommend Dilution: IHC 1:150

Reactivity: Human

Host: Mouse

Isotype: IgG1

Clonality: Monoclonal

**Immunogen:** Human recombinant protein fragment corresponding to amino acids 161-430 of human

NR0B1 (NP\_000466) produced in E.coli.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 51.5 kDa

**Gene Name:** nuclear receptor subfamily 0 group B member 1

Database Link: NP 000466 Entrez Gene 190 Human

**Background:** This gene encodes a protein that contains a DNA-binding domain. The encoded protein acts

as a dominant-negative regulator of transcription which is mediated by the retinoic acid receptor. This protein also functions as an anti-testis gene by acting antagonistically to Sry.

Mutations in this gene result in both X-linked congenital adrenal hypoplasia and

hypogonadotropic hypogonadism. [provided by RefSeq, Jul 2008]


Synonyms: AHC; AHCH; AHX; DAX-1; DAX1; DSS; GTD; HHG; NROB1; SRXY2

**Protein Families:** Druggable Genome, Nuclear Hormone Receptor, Transcription Factors






## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY NR0B1 ([RC204187], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-NR0B1. Positive lysates [LY424702] (100ug) and [LC424702] (20ug) can be purchased separately from OriGene.



Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human colon tissue using anti-NR0B1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris, pH8.5, 120°C for 3min, TA805006)