

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA805682

ALKBH1 Mouse Monoclonal Antibody [Clone ID: OTI1B8]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI1B8
Applications: IHC

Recommend Dilution: IHC 1:150

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 1-245 of human

ALKBH1 (NP_006011) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)


Gene Name: alkB homolog 1, histone H2A dioxygenase

Database Link: NP 006011 Entrez Gene 8846 Human

Synonyms: ABH; ABH1; alkB; ALKBH; hABH



Product images:



Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human colon tissue using anti-ALKBH1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris, pH8.5, 120°C for 3min, TA805682)