

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA806463

BCL2A1 Mouse Monoclonal Antibody [Clone ID: OTI3D10]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI3D10

Applications: WB

Recommend Dilution: WB 1:500
Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human BCL2A1 (NP_004040) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 20 kDa

Gene Name: BCL2 related protein A1

Database Link: NP 004040 Entrez Gene 597 Human

Background: This gene encodes a member of the BCL-2 protein family. The proteins of this family form

hetero- or homodimers and act as anti- and pro-apoptotic regulators that are involved in a wide variety of cellular activities such as embryonic development, homeostasis and

tumorigenesis. The protein encoded by this gene is able to reduce the release of proapoptotic cytochrome c from mitochondria and block caspase activation. This gene is a direct

transcription target of NF-kappa B in response to inflammatory mediators, and is upregulated by different extracellular signals, such as granulocyte-macrophage colony-stimulating factor (GM-CSF), CD40, phorbol ester and inflammatory cytokine TNF and IL-1, which suggests a cytoprotective function that is essential for lymphocyte activation as well as cell survival. Alternatively spliced transcript variants encoding different isoforms have been

found for this gene. [provided by RefSeq, Jul 2008]

Synonyms: ACC-1; ACC-2; ACC1; ACC2; BCL2L5; BFL1; GRS; HBPA1




BCL2A1 Mouse Monoclonal Antibody [Clone ID: OTI3D10] - TA806463

Protein Families: Druggable Genome

Protein Pathways: Metabolic pathways

Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY BCL2A1 ([RC201965], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-BCL2A1. Positive lysates [LY418253] (100ug) and [LC418253] (20ug) can be purchased separately from OriGene.