

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA806502**

## NARS2 Mouse Monoclonal Antibody [Clone ID: OTI5A9]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI5A9
Applications: IHC, WB

**Recommend Dilution:** WB 1:500, IHC 1:150

Reactivity: Human
Host: Mouse
Isotype: IgG2b

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human NARS2 (NP\_078954) produced in HEK293T

cell.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Predicted Protein Size:** 53.9 kDa

**Gene Name:** asparaginyl-tRNA synthetase 2, mitochondrial (putative)

Database Link: NP 078954 Entrez Gene 79731 Human


**Synonyms:** asnRS; DFNB94; SLM5

**Protein Pathways:** Aminoacyl-tRNA biosynthesis






## **Product images:**




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY NARS2 ([RC203201], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-NARS2 (1:500). Positive lysates [LY411139] (100ug) and [LC411139] (20ug) can be purchased separately from OriGene.




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human colon tissue using anti-NARS2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA806502) (1:150)






Immunohistochemical staining of paraffinembedded Carcinoma of Human lung tissue using anti-NARS2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA806502) (1:150)




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human endometrium tissue using anti-NARS2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA806502) (1:150)



Immunohistochemical staining of paraffinembedded Human lymphoma tissue using anti-NARS2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA806502) (1:150)





Immunohistochemical staining of paraffinembedded Human tonsil within the normal limits using anti-NARS2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA806502) (1:150)