

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA806565**

# SPC25 Mouse Monoclonal Antibody [Clone ID: OTI5E7]

### **Product data:**

Isotype:

**Product Type:** Primary Antibodies

Clone Name: OTI5E7
Applications: WB

Recommend Dilution: WB 1:500
Reactivity: Human
Host: Mouse

Clonality: Monoclonal

**Immunogen:** Full length human recombinant protein of human SPC25 (NP\_065726) produced in E.coli.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

lgG1

**Predicted Protein Size:** 26 kDa

**Gene Name:** SPC25, NDC80 kinetochore complex component

Database Link: NP 065726 Entrez Gene 57405 Human

**Background:** This gene encodes a protein that may be involved in kinetochore-microtubule interaction and

spindle checkpoint activity. [provided by RefSeq, Jul 2008]


Synonyms: AD024; hSpc25; SPBC25

**Protein Families:** Druggable Genome





# **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY SPC25 ([RC204623], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-SPC25. Positive lysates [LY412309] (100ug) and [LC412309] (20ug) can be purchased separately from OriGene.