

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA807146**

## Lymphocyte Activation Gene 3 (LAG3) Mouse Monoclonal Antibody [Clone ID: OTI10E7]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI10E7
Applications: FC, WB

Recommend Dilution: WB 1:2000

Reactivity: Human
Host: Mouse
Isotype: IgG2a

Clonality: Monoclonal

**Immunogen:** Human recombinant protein fragment corresponding to amino acids 66-332 of human

LAG3(NP\_002277) produced in E.coli.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Predicted Protein Size:** 54.8 kDa

**Gene Name:** lymphocyte activating 3

Database Link: NP 002277 Entrez Gene 3902 Human

**Background:** Lymphocyte-activation protein 3 belongs to Ig superfamily and contains 4 extracellular Ig-like

domains. The LAG3 gene contains 8 exons. The sequence data, exon/intron organization, and

chromosomal localization all indicate a close relationship of LAG3 to CD4. [provided by

RefSeq, Jul 2008]


Synonyms: CD223

**Protein Families:** Transmembrane






### **Product images:**




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY LAG3 ([RC220269], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-LAG3. Positive lysates [LY400828] (100ug) and [LC400828] (20ug) can be purchased separately from OriGene.



Western blot analysis of extracts (35ug) from 2 cell lines (Human PBMC treated with 10ug/ml PHA for 72h (Right)/untreated (Left)) lysates by using anti-LAG3 monoclonal antibody (1:150).





Flow cytometric analysis of living 293T cells transfected with LAG3 overexpression plasmid ([RC220269]), Red)/empty vector ([PS100001], Blue) using anti-LAG3 antibody (TA807146). Cells incubated with a non-specific antibody (Green) were used as isotype control. (1;100)