

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA807279

DR5 (TNFRSF10B) Mouse Monoclonal Antibody [Clone ID: OTI1B9]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI1B9

Applications: IHC, WB

Recommend Dilution: IHC 1:150

Reactivity: Human

Host: Mouse

Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 203-411 of human

TNFRSF10B(NP_671716) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 39.4 kDa

Gene Name: tumor necrosis factor receptor superfamily member 10b

Database Link: NP 671716 Entrez Gene 8795 Human

Background: The protein encoded by this gene is a member of the TNF-receptor superfamily, and contains

an intracellular death domain. This receptor can be activated by tumor necrosis factor-related apoptosis inducing ligand (TNFSF10/TRAIL/APO-2L), and transduces an apoptosis signal. Studies with FADD-deficient mice suggested that FADD, a death domain containing adaptor protein, is required for the apoptosis mediated by this protein. Two transcript variants encoding different isoforms and one non-coding transcript have been found for this gene.

[provided by RefSeg, Mar 2009]

Synonyms: CD262; DR5; KILLER; TRAIL-R2; TRAILR2; TRICK2; TRICK2A; TRICK2B; TRICKB; ZTNFR9

Protein Families: Druggable Genome, Transmembrane


Protein Pathways: Apoptosis, Cytokine-cytokine receptor interaction, Natural killer cell mediated cytotoxicity,

p53 signaling pathway






Product images:




Western blot analysis of extracts (10ug) from 4 different cell lines and 1 human tissue by using anti-TNFRSF10B monoclonal antibody at 1:200 dilution.




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human ovary tissue using anti-TNFRSF10B mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris, pH8.5, 120°C for 3min, TA807279)





Immunohistochemical staining of paraffinembedded Human pancreas tissue within the normal limits using anti-TNFRSF10B mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris, pH8.5, 120°C for 3min, TA807279)



Immunohistochemical staining of paraffinembedded Human prostate tissue within the normal limits using anti-TNFRSF10B mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris, pH8.5, 120°C for 3min, TA807279)