

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA807480

MRPL15 Mouse Monoclonal Antibody [Clone ID: OTI10B10]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI10B10
Applications: IHC, WB

Recommend Dilution: WB 1:500, IHC 1:500

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 28-296 of human

MRPL15(NP_054894) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 33.2 kDa

Gene Name: mitochondrial ribosomal protein L15

Database Link: NP 054894 Entrez Gene 29088 Human

Background: Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in

protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that

the latter contain a 5S rRNA. Among different species, the proteins comprising the

mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 39S subunit protein that belongs to the EcoL15 ribosomal protein family. A pseudogene corresponding to this


gene is found on chromosome 15q. [provided by RefSeq, Jul 2008]

Synonyms: HSPC145; L15mt; MRP-L7; MRP-L15; RPML7






Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY MRPL15 ([RC201471], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-MRPL15 (1:500). Positive lysates [LY415465] (100ug) and [LC415465] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 4 different cell lines by using anti-MRPL15 monoclonal antibody. (HeLa:human;293T:human;K562:human; MCF7:human) (1:500).






Immunohistochemical staining of paraffinembedded Human breast tissue within the normal limits using anti-MRPL15 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA807480) (1:500)




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human breast tissue tissue using anti-MRPL15 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA807480) (1:500)




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human ovary tissue using anti-MRPL15 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA807480) (1:500)






Immunohistochemical staining of paraffinembedded Human pancreas tissue within the normal limits using anti-MRPL15 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA807480) (1:500)




Immunohistochemical staining of paraffinembedded Carcinoma of Human thyroid tissue using anti-MRPL15 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA807480) (1:500)



Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human endometrium tissue using anti-MRPL15 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA807480) (1:500)





Immunohistochemical staining of paraffinembedded Human lymph node tissue within the normal limits using anti-MRPL15 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA807480) (1:500)