

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA807703

LIME (LIME1) Mouse Monoclonal Antibody [Clone ID: OTI5D3]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI5D3

Applications: IHC, WB

Recommend Dilution: WB 1:2000, IHC 1:150

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 30-295 of human

LIME1(NP_060276) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 31.1 kDa

Gene Name: Lck interacting transmembrane adaptor 1

Database Link: NP 060276 Entrez Gene 54923 Human

Background: LIME1 is a raft-associated transmembrane adaptor phosphoprotein that is preferentially

expressed in hemopoietic cells, particularly T cells (Brdickova et al., 2003 [PubMed 14610046];

Hur et al., 2003 [PubMed 14610044]). [supplied by OMIM, Mar 2008]


Synonyms: dJ583P15.4; LIME

Protein Families: Transmembrane






Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY LIME1 ([RC201909], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-LIME1 (1:2000). Positive lysates [LY413516] (100ug) and [LC413516] (20ug) can be purchased separately from OriGene.




Immunohistochemical staining of paraffinembedded Human prostate tissue within the normal limits using anti-LIME1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA807703) (1:150)





Immunohistochemical staining of paraffinembedded Carcinoma of Human prostate tissue using anti-LIME1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA807703) (1:150)



Immunohistochemical staining of paraffinembedded Human lymph node tissue within the normal limits using anti-LIME1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA807703) (1:150)



Immunohistochemical staining of paraffinembedded Human lymphoma tissue using anti-LIME1 mouse monoclonal antibody. (Heatinduced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA807703) (1:150)