

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA807792

NDUFB11 Mouse Monoclonal Antibody [Clone ID: OTI12H10]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI12H10

Applications: WB

Recommend Dilution: WB 1:2000
Reactivity: Human
Host: Mouse
Isotype: IgG2a

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 30-163 of human

NDUFB11(NP_061929) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 18.2 kDa

Gene Name: NADH:ubiquinone oxidoreductase subunit B11

Database Link: NP 061929 Entrez Gene 54539 Human

Background: NDUFB11 is a component of mitochondrial complex I. Complex I catalyzes the first step in the

electron transport chain, the transfer of 2 electrons from NADH to ubiquinone, coupled to the translocation of 4 protons across the membrane (Carroll et al., 2002 [PubMed 12381726]).

[supplied by OMIM, Feb 2009]


Synonyms: CI-ESSS; ESSS; Np15; NP17.3; P17.3

Protein Families: Transmembrane





Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY NDUFB11 ([RC213937], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-NDUFB11 (1:2000). Positive lysates [LY412795] (100ug) and [LC412795] (20ug) can be purchased separately from OriGene.