

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA807835

IGF1 Receptor (IGF1R) Mouse Monoclonal Antibody [Clone ID: OTI4C4]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI4C4
Applications: FC, WB

Recommend Dilution: WB 1:2000

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 741-935 of human

IGF1R(NP_000866) produced in Insect.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 81 kDa

Gene Name: insulin like growth factor 1 receptor

Database Link: NP 000866 Entrez Gene 3480 Human

Background: This receptor binds insulin-like growth factor with a high affinity. It has tyrosine kinase

activity. The insulin-like growth factor I receptor plays a critical role in transformation events. Cleavage of the precursor generates alpha and beta subunits. It is highly overexpressed in most malignant tissues where it functions as an anti-apoptotic agent by enhancing cell survival. Alternatively spliced transcript variants encoding distinct isoforms have been found

for this gene. [provided by RefSeq, May

Synonyms: CD221; IGFIR; IGFR; JTK13

Protein Families: Druggable Genome, Protein Kinase, Transmembrane

Protein Pathways: Adherens junction, Colorectal cancer, Endocytosis, Focal adhesion, Glioma, Long-term


depression, Melanoma, Oocyte meiosis, Pathways in cancer, Progesterone-mediated oocyte

maturation, Prostate cancer






Product images:





HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY IGF1R ([RC214928], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-IGF1R (1:2000). Positive lysates [LY400311] (100ug) and [LC400311] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 4 cell lines lysates by using anti-IGF1R monoclonal antibody (1:500).










Flow cytometric analysis of living HepG2 cells, using anti-IGF1R antibody (TA807835), (red), compared to an IgG isotype control, (green), and negative control (PBS), (blue) (1:25).



Flow cytometric analysis of living A549 cells, using anti-IGF1R antibody (TA807835), (red), compared to an IgG isotype control, (green), and negative control (PBS), (blue) (1:25).