

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA808318

RAB11FIP4 Mouse Monoclonal Antibody [Clone ID: OTI4C9]

Product data:


Product Type:	Primary Antibodies
Clone Name:	OTI4C9
Applications:	WB
Recommend Dilution:	WB 1:2000
Reactivity:	Human
Host:	Mouse
lsotype:	lgG2b
Clonality:	Monoclonal
Immunogen:	Human recombinant protein fragment corresponding to amino acids 91-177 of human RAB11FIP4(NP_116321) produced in E.coli.
Formulation:	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Predicted Protein Size:	71.7 kDa
Gene Name:	RAB11 family interacting protein 4
Database Link:	<u>NP_116321 Entrez Gene 84440 Human</u>
Background:	Proteins of the large Rab GTPase family (see RAB1A; MIM 179508) have regulatory roles in the formation, targeting, and fusion of intracellular transport vesicles. RAB11FIP4 is one of many proteins that interact with and regulate Rab GTPases (Hales et al., 2001 [PubMed 11495908]). [supplied by OMIM, Apr 2008]. ##Evidence-Data-START## Transcript exon combination :: AY169244.1, AJ314646.2 [ECO:0000332] RNAseq introns :: single sample supports all introns ERS025085, ERS025087 [ECO:0000348] ##Evidence-Data-END##
Synonyms:	RAB11-FIP4
Protein Pathways:	Endocytosis



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US



Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY RAB11FIP4 ([RC210877], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-RAB11FIP4 (1:2000). Positive lysates [LY409875] (100ug) and [LC409875] (20ug) can be purchased separately from OriGene.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US