

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA808440

TTF1 (NKX2-1) Mouse Monoclonal Antibody [Clone ID: OTI10F8]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI10F8

Applications: IF, IHC, WB

Recommend Dilution: IHC 1:150

Reactivity: Human

Host: Mouse

Isotype: IgG2b

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human NKX2-1(NP_003308) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 38.4 kDa

Gene Name: NK2 homeobox 1

Database Link: NP 003308 Entrez Gene 7080 Human

Background: This gene encodes a protein initially identified as a thyroid-specific transcription factor. The

encoded protein binds to the thyroglobulin promoter and regulates the expression of thyroid-specific genes but has also been shown to regulate the expression of genes involved in morphogenesis. Mutations and deletions in this gene are associated with benign hereditary chorea, choreoathetosis, congenital hypothyroidism, and neonatal respiratory distress, and may be associated with thyroid cancer. Multiple transcript variants encoding different isoforms have been found for this gene. This gene shares the symbol/alias 'TTF1' with another gene, transcription termination factor 1, which plays a role in ribosomal gene

transcription. [provided by RefSeq, Feb 2014]


Synonyms: BCH; BHC; EBP; NK-2; NKX2.1; NKX2A; NMTC1; T; TEBP; TITF1; TTF-1; TTF1

Protein Families: Druggable Genome, ES Cell Differentiation/IPS, Transcription Factors






Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY NKX2-1 ([RC217520], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-NKX2-1. Positive lysates [LY401141] (100ug) and [LC401141] (20ug) can be purchased separately from OriGene.




Immunohistochemical staining of paraffinembedded Human lung tissue within the normal limits using anti-NKX2-1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA808440) (1:150)






Immunohistochemical staining of paraffinembedded Carcinoma of Human lung tissue using anti-NKX2-1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA808440) (1:150)



Immunocytochemistry staining of NCI-H1975 cells using anti-NKX2-1 mouse monoclonal antibody (TA808440). The right is A549 cells as negative control.



Immunocytochemistry staining of TT cells using anti-NKX2-1 mouse monoclonal antibody (TA808440). The right is A549 cells as negative control.