

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA808577**

## **ERCC8 Mouse Monoclonal Antibody [Clone ID: OTI1G2]**

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI1G2
Applications: IHC, WB

**Recommend Dilution:** WB 1:500, IHC 1:150

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human ERCC8(NP\_001007235) produced in E.coli.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Predicted Protein Size:** 43.9 kDa

**Gene Name:** ERCC excision repair 8, CSA ubiquitin ligase complex subunit

Database Link: NP 001007235 Entrez Gene 1161 Human

Background: This gene encodes a WD repeat protein, which interacts with Cockayne syndrome type B

(CSB) protein and with p44 protein, a subunit of the RNA polymerase II transcription factor IIH. Mutations in this gene have been identified in patients with hereditary disease Cockayne syndrome (CS). CS cells are abnormally sensitive to ultraviolet radiation and are defective in the repair of transcriptionally active genes. Several transcript variants encoding different

isoforms have been found for this gene. [provided by RefSeq, Mar 2014]

Synonyms: CKN1; CSA; UVSS2


**Protein Families:** Druggable Genome, Transcription Factors

**Protein Pathways:** Nucleotide excision repair, Ubiquitin mediated proteolysis





## **Product images:**




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ERCC8 ([RC203124], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ERCC8 (1:500). Positive lysates [LY423459] (100ug) and [LC423459] (20ug) can be purchased separately from OriGene.



Immunohistochemical staining of paraffinembedded Human prostate tissue within the normal limits using anti-ERCC8 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA808577) (1:150)





Immunohistochemical staining of paraffinembedded Carcinoma of Human prostate tissue using anti-ERCC8 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA808577) (1:150)