

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA808629

GAS8 Mouse Monoclonal Antibody [Clone ID: OTI10G2]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI10G2

Applications: WB

Recommend Dilution: WB 1:2000

Reactivity: Human Host: Mouse

Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human GAS8 (NP_001472) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 56.2 kDa

Gene Name: growth arrest specific 8

Database Link: NP 001472 Entrez Gene 2622 Human

Background: This gene includes 11 exons spanning 25 kb and maps to a region of chromosome 16 that is

sometimes deleted in breast and prostrate cancer. The second intron contains an apparently

intronless gene, C16orf3, that is transcribed in the opposite orientation. This gene is a

putative tumor suppressor gene. Several transcript variants encoding different isoforms have


been found for this gene. [provided by RefSeq, Oct 2013]

Synonyms: CILD33; DRC4; GAS11





Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY GAS8 ([RC211150], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-GAS8 (1:2000). Positive lysates [LY419922] (100ug) and [LC419922] (20ug) can be purchased separately from OriGene.