

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA808680

HPCAL4 Mouse Monoclonal Antibody [Clone ID: OTI3C3]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI3C3

Applications: WB

Recommend Dilution: WB 1:2000

Reactivity: Human Host: Mouse

Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 2-191 of human

HPCAL4(NP_057341) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 22 kDa

Gene Name: hippocalcin like 4

Database Link: NP 057341 Entrez Gene 51440 Human

Background: The protein encoded by this gene is highly similar to human hippocalcin protein and

hippocalcin like-1 protein. It also has similarity to rat neural visinin-like Ca2+-binding proteintype 1 and 2 proteins. This encoded protein may be involved in the calcium-dependent

regulation of rhodopsin phosphorylation. The transcript of this gene has multiple

polyadenylation sites. Alternatively spliced transcript variants encoding different isoforms


have been found for this gene. [provided by RefSeq, Sep 2013]

Synonyms: HLP4





Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY HPCAL4 ([RC205253], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-HPCAL4 (1:2000). Positive lysates [LY414097] (100ug) and [LC414097] (20ug) can be purchased separately from OriGene.