

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA808796

PLCXD1 Mouse Monoclonal Antibody [Clone ID: OTI8C4]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI8C4

Applications: WB

Recommend Dilution: WB 1:2000

Reactivity: Human
Host: Mouse
Isotype: IgG2a

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human PLCXD1(NP_060860) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 36.5 kDa

Gene Name: phosphatidylinositol specific phospholipase C X domain containing 1

Database Link: NP 060860 Entrez Gene 55344 Human

Background: This gene is the most terminal protein-coding gene in the pseudoautosomal (PAR) region on

chromosomes X and Y. Alternate splicing results in multiple transcript variants. [provided by


RefSeq, Mar 2010]

Synonyms: LL0XNC01-136G2.1





Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY PLCXD1 ([RC202835], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-PLCXD1 (1:2000). Positive lysates [LY413089] (100ug) and [LC413089] (20ug) can be purchased separately from OriGene.