

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA808800

MAK Mouse Monoclonal Antibody [Clone ID: OTI5A7]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI5A7

Applications: WB

Reactivity: WB 1:500
Host: WB 0:500
WB 1:500
WB 1:500

Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 291-623 of human

MAK(NP_005897) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Gene Name: male germ cell associated kinase

Database Link: NP 005897 Entrez Gene 4117 Human

Background: The product of this gene is a serine/threonine protein kinase related to kinases involved in

cell cycle regulation. It is expressed almost exclusively in the testis, primarily in germ cells. Studies of the mouse and rat homologs have localized the kinase to the chromosomes during

meiosis in spermatogenesis, specifically to the synaptonemal complex that exists while

homologous chromosomes are paired. There is, however, a study of the mouse homolog that has identified high levels of expression in developing sensory epithelia so its function may be more generalized. Three transcript variants encoding different isoforms have been found for

this gene. [provided by RefSeq, Jul 2011]


Synonyms: dJ417M14.2; RP62

Protein Families: Druggable Genome, Protein Kinase





Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY MAK ([RC213532], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-MAK (1:500).