

### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA809199**

## P4HA3 Mouse Monoclonal Antibody [Clone ID: OTI9A2]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI9A2

Applications: IHC

Recommend Dilution: IHC 1:1000

Reactivity: Human Host: Mouse

Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombiant protein of human P4HA3 (NP\_878907) produced in E.coli.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Predicted Protein Size:** 59.3 kDa

Gene Name: prolyl 4-hydroxylase subunit alpha 3

Database Link: NP 878907 Entrez Gene 283208 Human

**Background:** This gene encodes a component of prolyl 4-hydroxylase, a key enzyme in collagen synthesis

composed of two identical alpha subunits and two beta subunits. The encoded protein is one of several different types of alpha subunits and provides the major part of the catalytic site of the active enzyme. In collagen and related proteins, prolyl 4-hydroxylase catalyzes the formation of 4-hydroxyproline that is essential to the proper three-dimensional folding of newly synthesized procollagen chains. Alternative splicing results in multiple transcript

variants. [provided by RefSeq, Jan 2014]

Synonyms: 2-oxoglutarate 4-dioxygenase (proline 4-hydroxylase); alpha III subunit; alpha polypeptide III;

C-P4H alpha III; collagen prolyl 4-hydroxylase alpha(III); procollagen-proline; prolyl 4-

hydroxylase


**Protein Families:** Druggable Genome

**Protein Pathways:** Arginine and proline metabolism, Metabolic pathways





# **Product images:**



Immunohistochemical staining of paraffinembedded Human skin tissue within the normal limits using anti-P4HA3 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA809199) (1:1000)