

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA809457

FOXO3 Mouse Monoclonal Antibody [Clone ID: OTI6H9]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI6H9

Applications: WB

Recommend Dilution: WB 1:2000

Reactivity: Human

Host: Mouse

Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 374-673 of human

FOXO3(NP_001446) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Gene Name: forkhead box O3

Database Link: NP 001446 Entrez Gene 2309 Human

Background: This gene belongs to the forkhead family of transcription factors which are characterized by a

distinct forkhead domain. This gene likely functions as a trigger for apoptosis through

expression of genes necessary for cell death. Translocation of this gene with the MLL gene is associated with secondary acute leukemia. Alternatively spliced transcript variants encoding

the same protein have been observed. [provided by RefSeq, Jul 2008]

Synonyms: AF6q21; FKHRL1; FKHRL1P2; FOXO2; FOXO3A

Protein Families: Druggable Genome, Transcription Factors


Protein Pathways: Chemokine signaling pathway, Endometrial cancer, Neurotrophin signaling pathway, Non-

small cell lung cancer





Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY FOXO3 ([RC202894], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-FOXO3 (1:500).