

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **Product datasheet for TA809528**

## WDR92 Mouse Monoclonal Antibody [Clone ID: OTI5H5]

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: OTI5H5

Applications: WB

Recommend Dilution: WB 1:2000

Reactivity: Human Host: Mouse

Isotype: IgG2b

Clonality: Monoclonal

Immunogen: Full length human recombiant protein of human WDR92 (NP\_612467) produced in HEK293T

cell.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Predicted Protein Size:** 39.6 kDa

**Gene Name:** WD repeat domain 92

Database Link: NP 612467 Entrez Gene 116143 Human

**Background:** This gene encodes a protein with two WD40 repeat domains thought to be involved in an

apoptosis via activation of caspase-3. Multiple transcript variants encoding different isoforms

have been found for this gene. [provided by RefSeq, Feb 2012]


Synonyms: FLJ31741

**Protein Families:** Druggable Genome





## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY WDR92 ([RC208985], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-WDR92 (1:2000). Positive lysates [LY408608] (100ug) and [LC408608] (20ug) can be purchased separately from OriGene.