

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA810001

N myc interactor (NMI) Mouse Monoclonal Antibody [Clone ID: OTI3C6]

Product data:

Applications:

Product Type: Primary Antibodies

WB

Clone Name: OTI3C6

Recommend Dilution: WB 1:2000

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 1-223 of human

NMI(NP_004679) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 34.9 kDa

Gene Name: N-myc and STAT interactor

Database Link: NP 004679 Entrez Gene 9111 Human

Background: NMYC interactor (NMI) encodes a protein that interacts with NMYC and CMYC (two members

of the oncogene Myc family), and other transcription factors containing a Zip, HLH, or HLH-Zip

motif. The NMI protein also interacts with all STATs except STAT2 and augments STAT-

mediated transcription in response to cytokines IL2 and IFN-gamma. The NMI mRNA has low

expression levels in all human fetal and adult tissues tested except brain and has high

expression in cancer cell line-myeloid leukemias. [provided by RefSeq, Jul 2008]


Synonyms: N-myc (and STAT) interactor; N-myc-interactor; N-myc and STAT interactor; N-myc interactor

Protein Families: Druggable Genome, Transcription Factors





Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY NMI ([RC201543], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-NMI (1:2000). Positive lysates [LY417813] (100ug) and [LC417813] (20ug) can be purchased separately from OriGene.