

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA810101

SUGT1 Mouse Monoclonal Antibody [Clone ID: OTI9A1]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI9A1 **Applications:** WB

Recommend Dilution: WB 1:500 Reactivity: Human Host: Mouse Isotype: IgG2a

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human SUGT1(NP_006695) produced in E.coli.

PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide. Formulation:

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 37.6 kDa

Gene Name: SGT1 homolog, MIS12 kinetochore complex assembly cochaperone

Database Link: NP 006695 Entrez Gene 10910 Human

Background: This gene is homologous to the yeast gene SGT1, which encodes a protein involved in

> kinetochore function and required for the G1/S and G2/M transitions. Complementation studies suggest that the human protein has similar functions. Two transcript variants

> encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2008]


Synonyms: SGT1

Protein Pathways: NOD-like receptor signaling pathway





Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY SUGT1 ([RC201676], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-SUGT1 (1:500). Positive lysates [LY416475] (100ug) and [LC416475] (20ug) can be purchased separately from OriGene.