

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# Product datasheet for TA810188

## PSMD9 Mouse Monoclonal Antibody [Clone ID: OTI5H2]

### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTI5H2
Applications:	WB
<b>Recommend Dilution:</b>	WB 1:500~2000
Reactivity:	Human
Host:	Mouse
lsotype:	lgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human PSMD9 (NP_002804) produced in HEK293T cell.
Formulation:	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Predicted Protein Size:	24.5 kDa
Gene Name:	proteasome 26S subunit, non-ATPase 9
Database Link:	<u>NP_002804 Entrez Gene 5715 Human</u>
Background:	The 26S proteasome is a multicatalytic proteinase complex with a highly ordered structure composed of 2 complexes, a 20S core and a 19S regulator. The 20S core is composed of 4 rings of 28 non-identical subunits; 2 rings are composed of 7 alpha subunits and 2 rings are composed of 7 beta subunits. The 19S regulator is composed of a base, which contains 6 ATPase subunits and 2 non-ATPase subunits, and a lid, which contains up to 10 non-ATPase subunits. Proteasomes are distributed throughout eukaryotic cells at a high concentration and cleave peptides in an ATP/ubiquitin-dependent process in a non-lysosomal pathway. An essential function of a modified proteasome, the immunoproteasome, is the processing of class I MHC peptides. This gene encodes a non-ATPase subunit of the 19S regulator. Three transcript variants encoding two different isoforms have been found for this gene. [provided]



View online »

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

by RefSeq, May 2012]

A549

#### Synonyms:


#### p27; Rpn4

170 \_\_\_\_\_ 130 \_\_\_\_ 100 \_\_\_\_ 70 \_\_\_\_ 55 \_\_\_

40 -

35 —


### **Product images:**



HELA HepG2 NCI-H460 Jurkat

HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY PSMD9 ([RC200903], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-PSMD9 (1:2000). Positive lysates [LY419093] (100ug) and [LC419093] (20ug) can be purchased separately from OriGene.

Western blot analysis of extracts (35ug) from 5 different cell lines by using anti-PSMD9 monoclonal antibody (1:500).



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US