

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA810205

DAPP1 Mouse Monoclonal Antibody [Clone ID: OTI2B5]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI2B5

Applications: WB

Recommend Dilution: WB 1:2000

Reactivity: Human Host: Mouse

Isotype: IgG2a

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human DAPP1 (NP_055210) produced in HEK293T

cell.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 32.6 kDa

Gene Name: dual adaptor of phosphotyrosine and 3-phosphoinositides 1

Database Link: NP 055210 Entrez Gene 27071 Human


Synonyms: BAM32

Protein Families: Druggable Genome, Phosphatase
Protein Pathways: B cell receptor signaling pathway





Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY DAPP1 ([RC215474], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-DAPP1 (1:2000). Positive lysates [LY402326] (100ug) and [LC402326] (20ug) can be purchased separately from OriGene.